

# Monoclonal Antibody to Human Interleukin-23 (p19)

BIOTINYLATED

<b>Antibody:</b>	MT155
<b>Product code:</b>	3457-6-250
<b>Quantity:</b>	250 µg
<b>Immunogen:</b>	Recombinant human IL-23
<b>Isotype:</b>	Mouse IgG2b
<b>Specificity:</b>	Native and recombinant human IL-23. The antibody recognises the p19 subunit of the IL-23 heterodimer (p40+p19). The antibody cross-reacts with native p19 from cynomolgus macaques. Please visit <a href="http://www.mabtech.com">www.mabtech.com</a> for information on reactivity with other monkey species.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Biotinylation:</b>	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of human IL-23 heterodimer (p40+p19) in solution e.g. cell culture supernatants using ELISA and enumeration of IL-23 producing cells using ELISpot. MAb MT155 is recommended as detection mAb in combination with coating mAbs MT86/221 (product code 3450-3).
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols.

# MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
[mabtech@mabtech.com](mailto:mabtech@mabtech.com)

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
[mabtech.usa@mabtech.com](mailto:mabtech.usa@mabtech.com)

Updated on, 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)