

Monoclonal antibody to Human EBI3

Antibody:	MT23A14
Product code:	3459-3-250
Quantity:	250 µg
Immunogen:	Recombinant human IL-27
Isotype:	Mouse IgG1
Specificity:	Native and recombinant human EBI3. MT23A14 binds to the EBI3 subunit of IL-27, IL-35, and IL-39.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of EBI3 in solution e.g. cell culture supernatants using ELISA and for enumeration of EBI3 producing cells using ELISpot. MT23A14 is recommended as coating mAb in combination with detection mAb MT140 (product code 3458-6).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com