

Monoclonal antibodies to Human Perforin

Antibody:	Pf-80/164 The reagent consists of a combination of two different monoclonal antibodies specific for Perforin.
Product code:	3465-3-1000
Quantity:	2000 µg
Immunogen:	Native human Perforin
Isotype:	Mouse IgG1
Specificity:	Native human Perforin. Cross-reacts with native Perforin from rhesus and cynomolgus macaques. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of Perforin in solution e.g. cell culture supernatants using ELISA and for enumeration of Perforin producing cells using ELISpot. Pf-80/164 is recommended as coating mAb in combination with detection mAb Pf-344 (product code 3465-6).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-04-07



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com