

Monoclonal Antibody to Human Perforin

BIOTINYLATED

Antibody:	Pf-344
Product code:	3465-6-1000
Quantity:	1000 µg
Immunogen:	Native human Perforin
Isotype:	Mouse IgG1
Specificity:	Native human Perforin. Cross-reacts with native Perforin from rhesus and cynomolgus macaques. Please visit www.mabtech.com for information on reactivity with other monkey species. The antibody also cross-reacts with Perforin from pig.
Concentration:	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of Perforin in solution e.g. cell culture supernatants using ELISA and for enumeration of Perforin producing cells using ELISpot. Pf-344 is recommended as detection mAb in combination with coating mAbs Pf-80/164 (product code 3465-3). Pf-344 is also suitable for immunocytochemistry, flow cytometry and Western Blot.
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Updated on, 2020-03-26

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com