

Monoclonal antibody to Human Interleukin-13

BIOTINYLATED

Antibody:	IL13-3
Product code:	3470-6-1000
Quantity:	1000 µg
Immunogen:	Recombinant human IL-13
Isotype:	Rat IgG1
Specificity:	Native and recombinant human IL-13. Cross-reacts with native IL-13 from rhesus and cynomolgus macaques and common marmosets. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Applications:	For quantification of IL-13 in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-13 producing cells using ELISpot. IL13-3 is recommended as detection mAb in combination with coating mAb IL13-I (product code 3470-3).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Intertek



Intertek

MABTECH AB
Box 1233
SE-131 28 Nacka Strand
Sweden
Tel +46 8 716 27 00
Fax +46 8 716 27 01
mabtech@mabtech.com

MABTECH Inc.
M.E.B. Suite 220, 3814 West Street
OH 45227 Cincinnati
USA
Tel: +1 513 871 4500
Fax: +1 513 871 7353
mabtech.usa@mabtech.com
Toll Free: +1 866 ELI-SPOT

MABTECH Australia Pty Ltd
Australia
Tel: +61 3 9470 4704
Fax: +61 3 8678 3216
mabtech.au@mabtech.com

MABTECH AB
Büro Deutschland
Germany
Tel: +49 40 4135 7935
Fax: +49 40 4135 7945
mabtech.de@mabtech.com

MABTECH AB
Bureau de liaison
France
Tel: +33 4 92 38 80 70
Fax: +33 4 92 38 80 71
mabtech.fr@mabtech.com