

# Monoclonal Antibody to Human Interleukin-5

BIOTINYLATED

<b>Antibody:</b>	5A10
<b>Product code:</b>	3490-6-250
<b>Quantity:</b>	250 µg
<b>Immunogen:</b>	Recombinant human IL-5
<b>Isotype:</b>	Rat IgG2a
<b>Specificity:</b>	Native and recombinant human IL-5. Cross-reacts with native IL-5 from rhesus and cynomolgus macaques. Please visit <a href="http://www.mabtech.com">www.mabtech.com</a> for information on reactivity with other monkey species.
<b>Concentration:</b>	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Biotinylation:</b>	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of IL-5 in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-5 producing cells using ELISpot. 5A10 is recommended as detection mAb in combination with coating mAb TRFK5 (product code 3490-3).
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols.

# MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tél: +46 8 716 27 00  
[mabtech@mabtech.com](mailto:mabtech@mabtech.com)

Updated on, 2020-03-26

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
[mabtech.usa@mabtech.com](mailto:mabtech.usa@mabtech.com)



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)