

to Human TNF- α

Antibody:	
Product code:	3510-3-1000
Quantity:	1000 μ g
Immunogen:	
Isotype:	
Specificity:	Native and recombinant human TNF- α . Cross-reacts with native TNF- α from rhesus and cynomolgus macaques. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of TNF- α in solution e.g. cell culture supernatants using ELISA and for enumeration of TNF- α producing cells using ELISpot. TNF3/4 are recommended as coating mAbs in combination with detection mAb TNF5 (product code 3510-6).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

Updated on 2020-03-26

www.mabtech.com