

Monoclonal Antibody to Monkey/Human TNF- α

BIOTINYLATED

Antibody:	MT15B15
Product code:	3512M-6-250
Quantity:	250 μ g
Immunogen:	Recombinant human TNF- α
Isotype:	Mouse IgG1
Specificity:	Native and recombinant human TNF- α . Cross-reacts with native TNF- α from rhesus and cynomolgus macaques. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of TNF- α in solution e.g. cell culture supernatants using ELISA and for enumeration of TNF- α producing cells using ELISpot. MT15B15 is recommended as detection mAb in combination with coating mAb MT21A8 (product code 3512M-3). The antibody is also suitable for Western blot and Flow cytometry.
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Updated on, 2020-03-26

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com