

Monoclonal Antibody to Human Interleukin-17A

Antibody:	MT44.6
Product code:	3520-3-1000
Quantity:	1000 µg
Immunogen:	Recombinant human IL-17A
Isotype:	Mouse IgG1
Specificity:	Native and recombinant human IL-17A. Cross-reacts with native IL-17 from African Green monkeys. Please visit www.mabtech.com for information on reactivity with other monkey species.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of IL-17A in solution e.g. cell culture supernatants using ELISA and for enumeration of IL-17A producing cells using ELISpot. MT44.6 is recommended as coating mAb in combination with detection mAb MT504 (product code 3520-6).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com