

# Monoclonal antibody to Human Latent TGF- $\beta$ 1

<b>Antibody:</b>	MT593
<b>Product code:</b>	3550-3-250
<b>Quantity:</b>	250 $\mu$ g
<b>Immunogen:</b>	Recombinant latent TGF- $\beta$ 1
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Native and recombinant latent Transforming Growth Factor- $\beta$ 1 (TGF- $\beta$ 1). The antibody binds to the human Latency Associated Protein (LAP), which is a part of the latent TGF- $\beta$ 1 complex. The antibody does not bind to human latent TGF- $\beta$ 2 or - $\beta$ 3. Cross-reacts with native latent TGF- $\beta$ 1 from rhesus and cynomolgus macaques.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 $\mu$ m) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of latent TGF- $\beta$ 1 in solution e.g. cell culture supernatants using ELISA. MT593 is recommended as coating mAb in combination with detection mAb MT517-biotin (product code 3550-6).
<b>Protocols:</b>	Please contact Mabtech for ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)