

# Monoclonal antibodies to Human Interleukin-8 (CXCL8)

<b>Antibody:</b>	MT8H6/8F19 The reagent consists of a combination of two different monoclonal antibodies specific for human IL-8.
<b>Product code:</b>	3560-1-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Recombinant human IL-8
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Native and recombinant human IL-8. Reacts with the 72 amino acid variant and 79 amino acid full-length human IL-8. Cross-reacts with IL-8 from cynomolgus macaque, IL-8 from cow and IL-8 from dog. Please visit <a href="http://www.mabtech.com">www.mabtech.com</a> for information on reactivity with other monkey species.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For enumeration of IL-8 producing cells using ELISpot. MT8H6/8F19 are recommended as coating mAbs in combination with detection mAb 26E5 (product code 3560-14).
<b>Protocols:</b>	Please contact Mabtech for ELISpot protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
[mabtech@mabtech.com](mailto:mabtech@mabtech.com)

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
[mabtech.usa@mabtech.com](mailto:mabtech.usa@mabtech.com)

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)