

Monoclonal Antibody to Human Interleukin-29 (IFN- λ 1)

BIOTINYLATED

Antibody:	MT6G4
Product code:	3570-6-250
Quantity:	250 μ g
Immunogen:	Recombinant human IL-29
Isotype:	Mouse IgG1
Specificity:	Native and recombinant human IL-29. Cross-reacts with human IL-28A and IL-28B.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 μ m) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of human IL-29 in solution, e.g. cell culture supernatants, using ELISA. MT6G4 is recommended as detection mAb in combination with IL-29 specific coating mAb MT2H7 (product code 3570-3).
Protocols:	Please contact Mabtech for ELISA protocols.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tél: +46 8 716 27 00
mabtech@mabtech.com

Updated on, 2020-03-26

Mabtech, Inc.
USA
Tél: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com