

Monoclonal antibody to Human CD43

Antibody:	B1B6
Product code:	3602-1-250
Quantity:	250 µg
Immunogen:	Human T cell leukemia
Isotype:	Mouse IgG1
Specificity:	Human CD43
Concentration:	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS. Contains no preservations.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	Studies of CD43 expression on cells and tissue sections by flow cytometry and immunostaining.
Protocols:	For futher information, please contact Mabtech.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com