

Monoclonal antibody to Human CD3

Antibody:	CD3-2
Product code:	3605-2
Quantity:	100 tests
Immunogen:	Stimulated human leucocytes
Isotype:	Mouse IgG2a
Specificity:	Recognizes the ϵ chain of human CD3
Contents:	Ready-to-use solution of conjugated antibody in sterile filtered (0.2 μ m) PBS with 0.2% BSA and 0.09% sodium azide. 5 μ l is recommended for staining of 1 million cells in a total volume of 50 μ l.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Conjugation:	The CD3-2 has been conjugated to fluorescein isothiocyanate (FITC).
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing.
Applications:	This mAb is useful for studies of CD3 expression in cells and tissue sections. For flow cytometry staining of paraformaldehyde fixed cells, fixation and permeabilisation are recommended after staining.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.



Intertek



Intertek

Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the above standards.

MABTECH AB
Box 1233
SE-131 28 Nacka Strand
Sweden
Tel +46 8 716 27 00
Fax +46 8 716 27 01
mabtech@mabtech.com

MABTECH Inc.
M.E.B. Suite 220, 3814 West Street
OH 45227 Cincinnati
USA
Tel: +1 513 871 4500
Fax: +1 513 871 7353
mabtech.usa@mabtech.com
Toll Free: +1 866 ELI-SPOT

MABTECH Australia Pty Ltd
Australia
Tel: +61 3 9470 4704
Fax: +61 3 8678 3216
mabtech.au@mabtech.com

MABTECH AB
Büro Deutschland
Germany
Tel: +49 40 4135 7935
Fax: +49 40 4135 7945
mabtech.de@mabtech.com

MABTECH AB
Bureau de liaison
France
Tel: +33 4 92 38 80 70
Fax: +33 4 92 38 80 71
mabtech.fr@mabtech.com