

Monoclonal antibody to Human CD8

Antibody:	CD8-1
Product code:	3607-1-1000
Quantity:	1000 µg
Immunogen:	Stimulated human leukocytes
Isotype:	Mouse IgG1
Specificity:	Recognises human CD8
Concentration:	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS. Contains no preservatives.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	Detection of CD8 on cells by e.g. immunofluorescent staining of cells or tissue sections. For flow cytometry staining of paraformaldehyde fixed cells, fixation is recommended after staining.
Protocols:	For further information, please contact Mabtech.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-04-07



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com