

Product Datasheet

ELISA substrate: pNPP for ALP

Product code:	3652-P10
Contents:	Contains pNPP, p-Nitrophenyl Phosphate, in diethanolamine buffer.
Applications:	Ready-to-use substrate for ELISA assays with alkaline phosphatase.
Instructions:	Add the substrate after incubation with alkaline phosphatase and subsequent washing. Incubate for typically 30-90 minutes at room temperature. The substrate produces a yellow product. Read the absorbance at 405-410 nm.
Storage:	Shipped at ambient temperature. Store at 4-8°C upon receipt.
Size:	120 ml



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2022-05-31



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com