

# Monoclonal antibody to Human Apolipoprotein A1

<b>Antibody:</b>	HDL 110
<b>Product code:</b>	3710-3-250
<b>Quantity:</b>	250 µg
<b>Immunogen:</b>	Human HDL
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Human apolipoprotein A1 (apoA1), purified or in the form of HDL. Cross-reacts with apoA1 from rhesus and cynomolgus macaques. The presence of low amounts of detergent is necessary for binding. Tween 20, Triton X100 or NP40 can be used at a concentration of 0.01-0.5%.
<b>Concentration:</b>	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of apoA1 in serum/ plasma samples and cell culture supernatants using ELISA. HDL 110 is recommended as coating mAb in combination with detection mAb HDL 44 (product code 3710-6). The antibody is also suitable for immunoprecipitation and Western blot.
<b>Protocols:</b>	Please contact Mabtech for ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)