

Monoclonal antibody to Human Apolipoprotein H

Antibody:	H219
Product code:	3711-3-1000
Quantity:	1000 µg
Immunogen:	Purified human apoH
Isotype:	Mouse IgG1
Specificity:	Human apolipoprotein H (apoH). Cross-reacts with apoH from rhesus and cynomolgus macaques.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of apoH in solution e.g. serum/plasma samples and cell culture supernatants using ELISA. H219 is recommended as coating mAb in combination with detection mAb H464 (product code 3711-6). The antibody is also suitable for Western blot.
Protocols:	Please contact Mabtech for ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com