

# Monoclonal Antibody to Human Apolipoprotein H

BIOTINYLATED

<b>Antibody:</b>	H464
<b>Product code:</b>	3711-6-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Purified human apoH
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Human apolipoprotein H (apoH). Cross-reacts with apoH from rhesus and cynomolgus macaques.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Biotinylation:</b>	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of apoH in solution e.g. serum/plasma samples and cell culture supernatants using ELISA. H464 is recommended as detection mAb in combination with coating mAb H219 (product code 3711-3). The antibody is also suitable for Western blot.
<b>Protocols:</b>	Please contact Mabtech for ELISA protocols.

# MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tél: +46 8 716 27 00  
mabtech@mabtech.com

Updated on, 2020-03-26

Mabtech, Inc.  
USA  
Tél: +1 513 871-4500  
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)