

Monoclonal antibody to Human Apolipoprotein J (Clusterin)

| | |
|-----------------------|--|
| Antibody: | J29 |
| Product code: | 3713-3-250 |
| Quantity: | 250 µg |
| Immunogen: | Purified human apoJ |
| Isotype: | Mouse IgG2b |
| Specificity: | Human apolipoprotein J (apoJ), Clusterin |
| Concentration: | Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide. |
| Purification: | Purified from <i>in vitro</i> cultures by protein G affinity chromatography. |
| Storage: | Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use. |
| Applications: | For quantification of apoJ in solution e.g. serum/plasma samples and cell culture supernatants using ELISA. J29 is recommended as coating mAb in combination with detection mAb J84 (product code 3713-6). The antibody is also suitable for Western blot. |
| Protocols: | Please contact Mabtech for ELISA protocols. |



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com