

Monoclonal Antibody to Human Apolipoprotein B

BIOTINYLATED

Antibody:	LDL 11
Product code:	3715-6-1000
Quantity:	1000 µg
Immunogen:	Human LDL
Isotype:	Mouse IgG1
Specificity:	Native human apoB in the form of LDL. Cross-reacts with LDL from rhesus and cynomolgus macaques. The presence of low amounts of detergent is necessary for binding. Tween 20, Triton X100 or NP40 can be used at a concentration of 0.01-0.5%.
Concentration:	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of apoB in serum/plasma samples and cell culture supernatants using ELISA. LDL 11 is recommended as detection mAb in combination with coating mAb LDL 20/17 (product code 3715-3). The antibody is also suitable for immunoprecipitation and Western blot.
Protocols:	Please contact Mabtech for ELISA protocols.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tél: +46 8 716 27 00
mabtech@mabtech.com

Updated on, 2020-03-26

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com