

# Monoclonal Antibody to Human Apolipoprotein B

<b>Antibody:</b>	LDL 11
<b>Product code:</b>	3715-7
<b>Quantity:</b>	250 µg
<b>Immunogen:</b>	Human LDL
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Native human LDL
<b>Contents:</b>	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Concentration:</b>	The LDL 11 has been conjugated to fluorescein isothiocyanate (FITC).
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing
<b>Applications:</b>	For intracellular staining of LDL using immunofluorescence. The simultaneous addition of fixative and a non-ionic detergent is recommended.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on 2021-10-05



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)