

Monoclonal antibody to Human Galectin-3

Antibody:	Gal-3
Product code:	3720-3-1000
Quantity:	1000 µg
Immunogen:	Mouse galectin-3
Isotype:	Rat IgG2a
Specificity:	Gal-3 detects galectin-3 in tissues and in cell suspensions.
Concentration:	Supplied at 1 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	Studies of antigen expression in cells and tissue sections using immunocytochemistry and immunoprecipitation.
Protocols:	For further information, please contact Mabtech.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-04-07



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com