

# Monoclonal antibody to Mouse Apolipoprotein A1

<b>Antibody:</b>	mHDL93
<b>Product code:</b>	3750-3-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Purified mouse apoA1
<b>Isotype:</b>	Rat IgG2a
<b>Specificity:</b>	Mouse apoA1
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of apoA1 in solution e.g. serum/plasma samples and cell culture supernatants using ELISA. mHDL93 is recommended as coating mAb in combination with biotinylated detection mAb mHDL36 (product code 3750-6). The antibody is also suitable for Western blot.
<b>Protocols:</b>	Please contact Mabtech for ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)