

Monoclonal antibody to Mouse Apolipoprotein A1

Antibody:	mHDL93
Product code:	3750-3-250
Quantity:	250 µg
Immunogen:	Purified mouse apoA1
Isotype:	Rat IgG2a
Specificity:	Mouse apoA1
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of apoA1 in solution e.g. serum/plasma samples and cell culture supernatants using ELISA. mHDL93 is recommended as coating mAb in combination with biotinylated detection mAb mHDL36 (product code 3750-6). The antibody is also suitable for Western blot.
Protocols:	Please contact Mabtech for ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com