

Monoclonal Antibody to Mouse Apolipoprotein E

BIOTINYLATED

Antibody:	mE29
Product code:	3752-6-1000
Quantity:	1000 µg
Immunogen:	Recombinant mouse apoE
Isotype:	Rat IgG2a
Specificity:	Mouse apoE
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For quantification of mouse apoE in serum/plasma samples and cell culture supernatants using ELISA. mE29 is recommended as detection mAb in combination with coating mAb mE1 (product code 3752-3). The antibody is also suitable for immunoprecipitation and Western blot.
Protocols:	Please contact Mabtech for ELISA protocols.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tél: +46 8 716 27 00
mabtech@mabtech.com

Updated on, 2020-03-26

Mabtech, Inc.
USA
Tél: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com