

# Monoclonal Antibody to Mouse IgE

BIOTINYLATED

<b>Antibody:</b>	MT44E
<b>Product code:</b>	3815-6-250
<b>Quantity:</b>	250 µg
<b>Immunogen:</b>	Mouse IgE
<b>Isotype:</b>	Rat IgG2b
<b>Specificity:</b>	Mouse IgE. The antibody is specific for mouse IgE allotypes IgE <sup>a</sup> (allotype a) and IgE <sup>b</sup> (allotype b).
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Biotinylation:</b>	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	<p>For quantification of IgE in serum, plasma and cell culture supernatants using ELISA and for enumeration of IgE producing cells using ELISpot. The antibody is also suitable for Western blot.</p> <p><b>IgE detection:</b> In ELISA and ELISpot, MT44E is recommended as detection mAb in combination with coating mAb MT56E (product code 3815-3)</p> <p><b>IgE<sup>a</sup> (allotype a) detection:</b> In ELISA and ELISpot, MT44E is recommended as detection mAb in combination with IgE<sup>a</sup> (allotype a) specific coating mAb MT18E (product code 3816-3).</p>
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols.

# MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on, 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)