

# Monoclonal antibody to Mouse IgE<sup>a</sup> (allotype a)

<b>Antibody:</b>	MT18E
<b>Product code:</b>	3816-3-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Mouse IgE
<b>Isotype:</b>	Rat IgG2a
<b>Specificity:</b>	Mouse IgE <sup>a</sup> (allotype a). The antibody is specific for mouse IgE <sup>a</sup> (allotype a) and does not bind IgE <sup>b</sup> (allotype b).
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Storage:</b>	Store product at 4-8°C. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For quantification of IgE <sup>a</sup> (allotype a) in serum, plasma and cell culture supernatants using ELISA and for enumeration of IgE <sup>a</sup> (allotype a) producing cells using ELISpot. MT18E is recommended as coating mAb in combination with biotinylated detection mAb MT44E (product code 3815-6).
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)