

# Monoclonal antibody to Mouse IgG2b

<b>Antibody:</b>	MG2b
<b>Product code:</b>	3823-3-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Mouse IgG
<b>Isotype:</b>	Rat IgG1
<b>Specificity:</b>	Mouse IgG2b. No cross-reactivity (<0.1%) with mouse IgG1, IgG2a, IgG2c (IgG2a analogue in e.g. C57BL/6 and NOD mice), IgG3, IgA, IgE, IgM or with IgG from rat, goat, rabbit, cow and human.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by ion-exchange chromatography.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For detection of mouse IgG2b in e.g. antigen-specific ELISA.
<b>Protocols:</b>	Please contact Mabtech for more information.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)