

# Monoclonal Antibody to Mouse IgG3

BIOTINYLATED

<b>Antibody:</b>	MG3
<b>Product code:</b>	3826-6-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Mouse IgG3
<b>Isotype:</b>	Rat IgG1
<b>Specificity:</b>	Mouse IgG3. No cross-reactivity (<0.1%) with mouse IgG1, IgG2a, IgG2b, IgG2c (IgG2a analogue in e.g. C57BL/6 and NOD mice), IgA, IgE, IgM or with IgG from rat, goat, rabbit, cow and human.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Biotinylation:</b>	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For detection of mouse IgG3 in e.g. antigen-specific ELISA.
<b>Protocols:</b>	Please contact Mabtech for more information.

# MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tél: +46 8 716 27 00  
mabtech@mabtech.com

Updated on, 2020-03-26

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)