

# Goat anti-mouse IgA

BIOTINYLATED

<b>Antibody:</b>	Antiserum is raised in goat using purified antigen as immunogen.
<b>Product code:</b>	3835-6-250
<b>Quantity:</b>	250 µg
<b>Immunogen:</b>	Mouse IgA
<b>Isotype:</b>	NA
<b>Specificity:</b>	The antibody is specific for mouse IgA H-chains. Negligible cross-reactivity with mouse IgG or IgM. May cross-react with immunoglobulins from other species.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Immunoaffinity chromatography
<b>Biotinylation:</b>	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
<b>Storage:</b>	Store product at 4-8°C. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	Enzyme immunoassays, Western blot techniques
<b>Protocols:</b>	Please contact Mabtech for ELISpot protocols.

# MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tél: +46 8 716 27 00  
mabtech@mabtech.com

Updated on, 2020-12-10

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)