

# Anti-human IgG mAb MT78; ALP conjugated

<b>Product code:</b>	3850-9A
<b>Antibody:</b>	MT78
<b>Enzyme:</b>	Highly purified enzyme from calf intestine with high activity especially suited for enzyme immunoassay purposes.
<b>Conjugation:</b>	The preparation has been thoroughly depleted of free ligand and free enzyme. The conjugate has high residual enzyme activity.
<b>Quantity:</b>	1 ml
<b>Applications:</b>	For quantification of IgG using e.g. ELISA. MAb MT78 is recommended as detection mAb in combination with coating mAb MT145 (product code 3850-3).
<b>Working dilution:</b>	Supplied in 0.1 M Tris buffer with 0.002% Kathon CG. Dilute 1/1000 as a general guideline. Dilution factor may vary with differing assay conditions.
<b>Storage:</b>	Shipped at ambient temperature. On arrival store product at 4-8°C.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tel: +46 8 716 27 00  
mabtech@mabtech.com

Updated on 2020-04-02

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)