

# Monoclonal Antibody to Human IgG1

BIOTINYLATED

<b>Antibody:</b>	MTG1218
<b>Product code:</b>	3851-14-1000
<b>Quantity:</b>	1000 µg
<b>Immunogen:</b>	Human IgG
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Human IgG1. No cross-reactivity (<0.1%) with human IgG2, IgG3, IgG4, IgA1, IgA2, IgE, IgM, IgD or with IgG from mouse and cow. Cross-reacts with IgG1 from rhesus and cynomolgus macaque. Please visit <a href="http://www.mabtech.com">www.mabtech.com</a> for information on reactivity with other monkey species.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Biotinylation:</b>	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For detection of IgG1 in e.g., antigen-specific ELISA and ELISpot.
<b>Protocols:</b>	Please contact Mabtech for ELISpot and ELISA protocols.

# MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tél: +46 8 716 27 00  
[mabtech@mabtech.com](mailto:mabtech@mabtech.com)

Updated on, 2020-05-04

Mabtech, Inc.  
USA  
Tél: +1 513 871-4500  
[mabtech.usa@mabtech.com](mailto:mabtech.usa@mabtech.com)



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)