

# Monoclonal Antibody to Human IgG3

BIOTINYLATED

<b>Antibody:</b>	MTG34
<b>Product code:</b>	3853-6-250
<b>Quantity:</b>	250 µg
<b>Immunogen:</b>	Human IgG3
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	Human IgG3. No cross-reactivity (<0.1%) with human IgG1, IgG2, IgG4, IgA, IgE, IgM or with IgG from mouse, rat, goat, rabbit and cow. Low cross-reactivity (<1%) with IgG from rhesus and cynomolgus macaque. Please visit <a href="http://www.mabtech.com">www.mabtech.com</a> for information on reactivity with other monkey species.
<b>Concentration:</b>	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
<b>Purification:</b>	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
<b>Biotinylation:</b>	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
<b>Storage:</b>	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
<b>Applications:</b>	For detection of IgG3 in e.g. antigen-specific ELISA.
<b>Protocols:</b>	Please contact Mabtech for more information.

# MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)  
Sweden  
Tél: +46 8 716 27 00  
[mabtech@mabtech.com](mailto:mabtech@mabtech.com)

Updated on, 2020-03-26

Mabtech, Inc.  
USA  
Tel: +1 513 871-4500  
[mabtech.usa@mabtech.com](mailto:mabtech.usa@mabtech.com)



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

[www.mabtech.com](http://www.mabtech.com)