

Monoclonal Antibody to Mouse IgA

BIOTINYLATED

Antibody:	MT39A
Product code:	3865-6-1000
Quantity:	1000 µg
Immunogen:	Mouse immunoglobulins
Isotype:	Rat IgG1
Specificity:	Mouse IgA
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Biotinylation:	Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	For enumeration of IgA (total and antigen-specific) producing cells using ELISpot and for quantification of antigen-specific IgA using e.g. indirect ELISA. In ELISpot, MT39A is recommended as detection mAb in combination with coating mAb MT45A (product code 3865-3).
Protocols:	Please contact Mabtech for ELISpot and ELISA protocols.

MABTECH

Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tél: +46 8 716 27 00
mabtech@mabtech.com

Updated on, 2020-03-26

Mabtech, Inc.
USA
Tél: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com