

Monoclonal antibody to Mouse Ig kappa chain

Antibody:	MT24
Product code:	3875-3-1000
Quantity:	1000 µg
Immunogen:	Mouse immunoglobulins
Isotype:	Rat IgG1
Specificity:	Mouse Ig kappa light chain. No cross-reactivity (<0.1%) with Ig from rat, goat, rabbit, cow, sheep, horse or human.
Concentration:	Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide.
Purification:	Purified from <i>in vitro</i> cultures by protein G affinity chromatography.
Storage:	Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use.
Applications:	Antigen-specific ELISA and ELISpot
Protocols:	Please contact Mabtech for more information.



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tel: +46 8 716 27 00
mabtech@mabtech.com

Mabtech, Inc.
USA
Tel: +1 513 871-4500
mabtech.usa@mabtech.com

Updated on 2020-03-26



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com