

Monoclonal Antibody to Mouse Ig kappa chain

BIOTINYLATED

| | |
|-----------------------|---|
| Antibody: | MT24 |
| Product code: | 3875-6-250 |
| Quantity: | 250 µg |
| Immunogen: | Mouse immunoglobulins |
| Isotype: | Rat IgG1 |
| Specificity: | Mouse Ig kappa light chain. No cross-reactivity (<0.1%) with Ig from rat, goat, rabbit, cow, sheep, horse or human. |
| Concentration: | Supplied at 0.5 mg/ml in sterile filtered (0.2 µm) PBS with 0.02% sodium azide. |
| Purification: | Purified from <i>in vitro</i> cultures by protein G affinity chromatography. |
| Biotinylation: | Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin. |
| Storage: | Store product at 4-8°C or frozen at -20°C or below. Avoid repeated freezing/thawing. The expiry date indicates how long unopened products, stored according to instructions, are recommended for use. |
| Applications: | Antigen-specific ELISA and ELISpot |
| Protocols: | Please contact Mabtech for more information. |



Note; for research use only.

Mabtech shall not be liable for the use or handling of the product or for consequential, special, indirect or incidental damages therefrom.

Mabtech AB (Head Office)
Sweden
Tél: +46 8 716 27 00
mabtech@mabtech.com

Updated on, 2020-03-26

Mabtech, Inc.
USA
Tél: +1 513 871-4500
mabtech.usa@mabtech.com



Developed and manufactured by MABTECH AB, Sweden, whose quality management system complies with the standards ISO 9001:2015 & ISO 13485:2016.

www.mabtech.com